COOK
MEDICAL

CE 0123

**EN
3**

**Biodesign®
Surgisis® 8-Layer Tissue Graft
Surgisis® Hernia Graft**

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
Bruksanvisning

FP0036-02J

BIODESIGN® SURGISIS® 8-LAYER TISSUE GRAFT SURGISIS® HERNIA GRAFT

INTENDED USE

The Cook® Biodesign® Surgisis® Hernia Graft and the 8-Layer Tissue Graft are intended to be implanted to reinforce soft tissues where weakness exists. Indications for use include the repair of a hernia or body wall defect. The graft is supplied sterile and is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

SURGISIS® HERNIA GRAFT This symbol means the following:
Surgisis® Hernia Graft

SURGISIS® 8-LAYER TISSUE GRAFT This symbol means the following:
Surgisis® 8-Layer Tissue Graft

These grafts are intended for use by trained medical professionals.

CONTRAINDICATIONS

These grafts are derived from a porcine source and should not be used in patients with known sensitivity to porcine material.

PRECAUTIONS

- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- **Do not resterilize.** Discard all open and unused portions of the graft.
- The graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard graft if mishandling has caused possible damage or contamination, or if the graft is past its expiration date.
- Ensure that graft is rehydrated prior to cutting, suturing, stapling, tacking or loading of the graft laparoscopically.
- Ensure that all layers of the graft are secured when suturing, stapling, or tacking.
- Place graft in maximum possible contact with healthy, well-vascularized tissue to encourage cell ingrowth and tissue remodeling.
- Suturing, stapling, or tacking more than one graft together may decrease graft performance.
- No studies have been conducted to evaluate the reproductive impact of the clinical use of the graft.
- Extended rehydration or excessive handling could lead to partial delamination of superficial layers of the graft.
- This graft has been reported to be safe in clean and clean-contaminated hernia repair.¹ However, care should be exercised when it is placed in critically ill patients or in severely contaminated wounds.
- Care should be taken to avoid damage to the graft when loading laparoscopically. It is recommended to load through a 10 mm or larger port.
- When using with negative pressure vacuum-assisted devices, cover the graft with a hydrogel or non-adherent dressing to prevent dryness. Additionally, use negative pressure wound therapy systems according to manufacturer recommendations including prevention of air leaks and dressing changes.

POTENTIAL COMPLICATIONS

Possible adverse reactions with the use of any prosthesis may include, but are not limited to:

- infection
- inflammation
- adhesion
- fistula formation
- seroma formation
- hematoma
- bowel erosion
- recurrence of tissue defect
- premature degradation

Complications, such as delayed wound infection, premature degradation, hernia recurrence, bowel erosion, and the need for re-operation, should be reasonably expected in patients who are critically ill or who have severely contaminated abdomens.

STORAGE

This graft should be stored in a clean, dry location at room temperature.

STERILIZATION

This graft has been sterilized with ethylene oxide.

INSTRUCTIONS FOR USE

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

NOTE: Always handle the graft using aseptic technique, minimizing contact with latex gloves.

REQUIRED MATERIALS

- Sterile forceps
 - Rehydration fluid: room temperature, sterile saline or sterile lactated Ringer's solution.
1. Using aseptic technique, remove the sterile inner package from the outer packaging and place it in the sterile field.
 2. Open the sterile inner package carefully. Aseptically, remove the graft from the inner package.
 3. Rehydrate the graft in room temperature, sterile saline or sterile lactated Ringer's solution until the desired handling characteristics are achieved. For the Hernia Graft, rehydrate for less than one (1) minute. For the 8-Layer Tissue Graft, rehydrate for at least five (5) minutes prior to cutting, suturing, stapling, tacking or loading the device laparoscopically. If the graft is provided in a tray, it can be rehydrated in the tray.
 4. Prepare the hernia repair site using standard surgical techniques.

MESH PLACEMENT

- Preperitoneal/retrorectus (Rives and Stoppa technique) or intraperitoneal placement (underlay) is recommended for open procedures.
 - Intraperitoneal placement (underlay) is recommended for a laparoscopic approach.
5. Using aseptic technique, trim the graft to fit the site, providing an allowance for overlap. Fundamental surgical principles suggest recurrence can be minimized if the mesh overlaps the surrounding tissue by at least 4-5 cm in all directions.²

NOTE: If the graft is cut too small for the defect, excess tension may be placed on the suture line. This can result in recurrence of the original tissue defect or development of a defect in the adjacent tissues.

6. Using aseptic technique, transfer the graft to the surgical site and suture, staple, or tack into place, avoiding excess tension.

NOTE: Excessive handling may lead to partial delamination of a superficial layer. This can be managed by securing the detached portion to adjacent tissue with sutures or staples.

NOTE: Surgical experience indicates that suturing or stapling the graft with close tissue approximation produces better outcomes. Leave less than or equal to 3cm between sutures with bite depth of 1cm. Use permanent or long-term absorbable sutures or tacks according to surgeon preference.

NOTE: Interrupted sutures can provide additional security against recurrence of tissue defect in the event of suture failure.

NOTE: The Hernia Graft has been designed with an open matrix structure to optimize tissue incorporation and angiogenesis.

NOTE: It is recommended to minimize the bridging of defects to improve outcomes.

7. Complete the standard surgical procedure.
8. Discard any unused portions according to institutional guidelines for medical waste.
9. Place closed suction drains for 2-6 weeks. Remove when output is less than 20 mL/24 hours for at least two (2) consecutive days or until drain is dry.

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© COOK BIOTECH 2013
West Lafayette, IN 47906 U.S.A.
AUSTRALIA PTY. LTD.